

Scientific & Educational Software

Clone Manager File Viewer

The Clone Manager File Viewer for Windows Mobile Devices, version 1.0, can be used on pocket PC's or PDA's running Windows Mobile ver 5, Windows Mobile 2003, Windows Mobile 2002 or Windows CE 3.0. A stylus is required for option selection.

Installation Instructions

You can use the downloaded setup program (CMViewer.msi) to install the Clone Manager File Viewer on your desktop computer using the Windows Installer for subsequent download of this application to your mobile device using Active Sync. To start the Clone Manager File Viewer installation:

- Connect your Pocket PC or PDA to your desktop computer.
- Double-click on the downloaded CMViewer.msi file, saved to your computer hard disk.
- Follow the instructions on your screen.

The installation program will copy the program files to the specified installation directory on your desktop computer. When installation is complete, Active Sync will load the application into Program Files on your mobile device. If you are not able to use Active Sync, you can install the program to your PC to generate CAB files and then copy the appropriate CAB file to your mobile device (see instructions on next page).

To run the application, tap Start > Programs > CMViewer.

Molecule Files

The viewer can display Clone Manager molecule files (*.cm5 or *.cx5) or primer collection files (*.px5). You can use Active Sync to move copies of your files to your mobile device or simply copy files to your storage card. To configure Active Sync to transfer your files, start Active Sync on your desktop and click Tools, Options and set the checkbox to enable transfer of Files. Active Sync will set up a folder on your desktop where you can place the files you want to transfer during synchronization. After synchronization, the molecule files can be found in the My Documents folder on the mobile device.

User Documentation

A User Guide for the Clone Manager File Viewer, version 1.0 can be viewed or printed using the Adobe Acrobat reader. This document is available at the CMViewer download page at the Sci Ed Software web site. While you are using the Clone Manager File Viewer on your mobile device, you can tap Tools > Help to access application help topics.

Uninstall Instructions

To remove the application from your desktop computer, click Start, Settings, Control Panel and then Add/Remove Programs. Select CMViewer and click Remove button. If installed through Active Sync, the program will also be removed from the mobile device, if connected. To remove from your mobile device, tap Start>Settings>System>Remove Programs. Select Sci Ed Software CMViewer and tap remove.

Alternate Installation Instructions

If you are unable to use Active Sync, you can install the Clone Manager File Viewer for Mobile Devices by copying an appropriate CAB file to your Pocket PC or PDA.

To install:

- Double-click on the downloaded CMViewer.msi file, saved to your computer hard disk.
- Follow the instructions on your screen to install the program on your desktop computer.
- Use File Explorer to navigate to the installation directory and locate the CAB files. (The installation directory is usually Program Files\Sci Ed Software\CMViewer\Setup).
- Select the appropriate CAB file, using the recommendations below.
- Copy this file to your mobile device or the storage card used for this device.
- On the mobile device, use File Explorer to locate the CAB file you just copied.
- Click on this file to install the Clone Manager File Viewer on your device.

CAB Files:

CMViewer_PPC.ARMV4.CAB -- For Pocket PC2003 Strong Arm and XScale processors, and Windows Mobile 5 Intel XScale processors.

CMViewer_PPC.ARM.CAB – for Pocket PC 2000 and 2002 StrongArm and XScale processors.

CVIEWER_PPC.MIPS.CAB, CMViewer_PPC.SH3.CAB, CMViewer_PPC.WCE220X86.CAB – other processors, usually legacy mobile devices.