

Mobile Viewers

Installation Instructions

The mobile viewers are cross platform apps that run on a variety of devices including Android, Apple and Windows Store smart phones and tablets. Installation uses the appropriate platform's app, or play, store. More details are given below.

- Open the appropriate App, or Play, Store on your device
- Search for 'SciEd Clone Manager Viewer'
- Select the option to install

The installer will install the viewer app and create an icon to launch the app. A sample collection and some molecules will be installed.

Tutorials

Tutorials covering using the mobile viewer apps and sharing molecules with the Clone Manager desktop application are available in the tutorials section of our web site (https://scied.com/doc_tutorials11.htm)